



# SAVE II+ VENTILATOR

The SAVE II+ ventilator is a small, extremely durable, portable mechanical ventilator designed to provide life-saving mechanical ventilation in pre-hospital, aeromedical, and field hospitals. The SAVE II+ is not a full-featured ICU ventilator; however, it can extend response capabilities for surge requirements during a pandemic or mass-casualty situation for adult patients weighing at least 45kg.

## TACTICAL COMBAT CASUALTY CARE CONFIGURATION

Built for the battlefield, the SAVE II+ is durable, portable, and powered for up to 9.25 hours by battery or AC power. The SAVE II+ is designed for trauma medics where size, weight, and extreme ease of use are paramount. Simply select the patient's height, and the device dials in a preset tidal volume of 6mL/kg of the patient's ideal body weight. This can help reduce operator error and eliminates the guesswork associated with bagging in a high-stress environment or setting up an overly-complicated transport vent.

## IN ANY ENVIRONMENT

The ease of use, portability, and long battery life make the SAVE II+ well-suited to provide life-saving mechanical ventilation to the developing world, field hospitals, or other environments without ICU equipment. The design of the SAVE II+ enables safe and effective use with minimal training.

Preset buttons provide ARDSNet-compliant rate/volume based on the ideal body weight for a given patient.

Adjustable PIP (10–60) and PEEP (0–20).

Supports minute ventilation up to 12.5 LPM.

Accepts low-pressure supplemental oxygen to provide up to 100% FiO<sub>2</sub>.

Runs for up to 9.25 hours per full battery charge.

Detects blockages and disconnects.

Durable, lightweight, compact design.



visit [safeguardmedical.com](http://safeguardmedical.com)  
call (US): +1 (855) 428-6074 / (Int'l): +44 (0) 1568 613942  
email (US): [info@safeguardmedical.com](mailto:info@safeguardmedical.com) / (Int'l): [sales.admin@safeguardmedical.com](mailto:sales.admin@safeguardmedical.com)



**SPECIFICATIONS**

PN	70110H	NSN	6515-01-696-5218
Device dims.	6.5" x 6.25" x 2" (16.5 x 15.9 x 5.1cm)	Device weight	2.8lbs (1.3kg)
Case dims.	14" x 16" x 7" (35.6 x 40.6 x 17.8cm)	Case weight	9.9lbs (4.5kg)

**DESIRED SURGE PERFORMANCE CHARACTERISTICS<sup>1</sup>**

**SAVe II+ PERFORMANCE CHARACTERISTICS**

Dimensions	Small enough to carry	6.5" x 6.25" x 2" (16.5cm x 15.9cm x 5.1cm)
Weight	<22lbs (10kg)	2.8lbs (1.3kg)
Flow rate	Not specified	0–40 LPM
I:E ratio	Not specified	1:2
Respiratory rate	6–35 BPM	8–30 BPM
Tidal volume	250–750mL	200–800mL
Peak inspiratory pressure (cmH <sub>2</sub> O)	<30	10–60
Peak end expiratory pressure (cmH <sub>2</sub> O)	5–15	0–20
Battery duration	>4 hours	8.25–9.25 hours
FiO <sub>2</sub>	21%–100%	✓
Controls	<ul style="list-style-type: none"> <li>• Separate TV &amp; RR control</li> <li>• Monitor airway pressure</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> </ul>
Meets ARDSNet guidelines	Yes	✓
Price (USD)	<\$10,000	✓
Mode	Volume control verification	✓
Country	USA	✓
Alarms	<ul style="list-style-type: none"> <li>• Circuit disconnect</li> <li>• High airway pressure</li> <li>• Low airway pressure</li> <li>• Loss of electrical power</li> <li>• Loss of high-pressure gas</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> </ul>
Ease of use	Body-weight prediction based on height	✓

1. Surge Capacity Mechanical Ventilation, Respiratory Care, Jan 2008 Vol 53 No 1, Branson et al



visit [safeguardmedical.com](http://safeguardmedical.com)  
 call (US): +1 (855) 428-6074 / (Int'l): +44 (0) 1568 613942  
 email (US): [info@safeguardmedical.com](mailto:info@safeguardmedical.com) / (Int'l): [sales.admin@safeguardmedical.com](mailto:sales.admin@safeguardmedical.com)

