

Pelvic Splint



Overview

The Pelvic Splint consists of a central strap made of neoprene with two straps running through a central buckle. When aligned around the pelvis and a compressive force is applied, stabilisation of the pelvis can be achieved.

Key Benefits

- A neoprene band with no internal stitching, which aims to reduce the risk of pressure ulceration and tissue damage.
- Tightening straps are fastened to the neoprene using secure hook material.
- Device is one size fits all, and can be trimmed using scissors.
- The device can be applied and adjusted by a single operator.
- The device allows access to the abdomen and groin for urinary catheterization and for accessing the femoral vein.
- The device is single use (disposable).
- Product code: PR0400 (Standard) / PRO410 (Tactical)

Specifications



Dimensions: 1260mm x 190mm



Weight: Under 250g



Product shelf life: 5 years

For further information on the Pelvic Splint, visit: www.safeguardmedical.com