



NASAL AIRWAY

Airway management in a pre-hospital setting can be challenging, and is critical to patient survival. A nasopharyngeal airway (NPA) is a simple-to-insert medical device, used to maintain an open airway by responders at all levels. The unique flange serves as an external bevel indicator and is a securing feature to hold the device in place.

The NPA is constructed of medical-grade silicone which resists the effects of extreme temperatures, retains no memory from storage, and maintains its shape, ensuring an open path for air passage after insertion. The NPA is sized to 28 French, or size 7, to fit the small-to-medium-sized adult population for easy first responder program fielding, and a water-soluble jelly sachet is included to facilitate proper insertion.¹

FEATURES & BENEFITS

- **Flange & Rounded Bevel Tip:** Allows for gentle insertion, securing, and bevel indication
- **Silicone Material:** Will not kink, crimp, or lose shape, and avoids damage to the nasal passage while maximizing patient comfort
- **Lightweight & Low Cube:** Packaged for ease of storage
- **Water-soluble jelly sachet** included for insertion¹
- **Easy identification** with color-coded labeling

1. Atanelov Z, Aina T, Amin B, et al. Nasopharyngeal Airway. [Updated 2022 Sep 26]. In: StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2022 Jan-. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK513220/>

Devices included in this document may not be cleared or licensed for use or sale in your individual country. Contact your local distributor for availability. SGM-MKT-NAW-01-Rev-002 | © 2026 Safeguard US Operating, LLC | 02/20/2026

SPECIFICATIONS

- **Dimensions[†]** **10.8" x 4.5"**
(27.3cm x 11.4cm)
- **Weight** **0.75oz**
(21g)
- **Shelf life*** (NPA) **6 years**

PART NUMBERS

- **11-128**
Pre-Lubricated Nasal Airway (sterile)
Availability: US, Singapore
- **PDF129**
Nasal Airway (w/ lubricant sachet)
Availability: US, Saudi Arabia, Brazil
- **PDF228**
28Fr Nasal Airway (w/o lubricant)
Availability: US, UK, EU
- **PP0054** (Procedure Pack)
28Fr Nasal Airway (w/ lubricant)
Availability: UK, EU

[†]Packaged | ^{*}From date of manufacture