



TMT

TACTICAL MECHANICAL TOURNIQUET™

CoTCCC-RECOMMENDED

The CoTCCC-recommended (Committee on Tactical Combat Casualty Care) TMT—or Tactical Mechanical Tourniquet™—is a hemorrhage control device specifically designed for massive hemorrhage control of an extremity.

Using lessons learned from the battlefield, the TMT ensures ease of application and definitive hemorrhage control, and requires minimal training to instantly treat life-threatening hemorrhage of an extremity.

DESIGNED TO SAVE LIVES

INCREASED VISIBILITY

The high-visibility TMT ensures rapid identification and application in critical, high-stress situations, significantly improving response time when every second counts.

DUAL LOCKING MECHANISMS

Double-routed speed buckle and hook-and-loop backed strap prevents slippage and pressure loss.

ULTRA-DURABLE

A UV-protective additive to polymers means the TMT is durable enough to be stored outside of its packaging.

LOWER OCCLUSION PRESSURE

2-inch wide strap achieves occlusion at lower pressures, while quick-disconnect buckle prevents misrouting of tail and allows for trapped-limb or one-handed application.



DETAILS & SPECS

PN	PDF186
Packaged	11.4cm x 5.1cm x 6.2cm (4.5" x 2" x 2.5")
Weight	100g (3.5oz) nominal as packaged
Shelf life	7 years from date of manufacture
Accommodates limb circumferences	5.75"–38" (14.6cm–96.5cm)

