



CHITOGAUZE® XR PRO

Across the globe, ChitoGauze® XR PRO is a top choice in hemostatic gauze amongst first responders. This chitosan-impregnated product provides life-saving actions when seconds matter¹ and is derived from a naturally occurring, biocompatible polysaccharide that is designed to control moderate to severe hemorrhage. The hemostatic effect is achieved via positively-charged chitosan interacting with negatively-charged red blood cell membranes, resulting in a life-saving clot.

ChitoGauze XR PRO readily conforms to a wide range of wound profiles, allowing for easy wound packing and rapid control of blood loss. ChitoGauze XR PRO packs wounds like a traditional user-friendly gauze and features a radio-opaque strip.

In a world of many hemostatic agents, ChitoGauze XR PRO is the medic's choice for your medical kit or aid bag.

1. Te Grotenhuis, R., van Grunsven, P. M., Heutz, W. M., & Tan, E. C. (2016). Prehospital use of hemostatic dressings in emergency medical services in the Netherlands: A prospective study of 66 cases. *Injury*, 47(5), 1007-1011. <https://doi.org/10.1016/j.injury.2016.01.005>



REDUCED SHEDDING:

Chitosan-impregnated gauze specifically designed to limit particle shedding into the wound.

RADIO-OPAQUE STRIP: Easy to identify on radiographic imaging.

Z-FOLDED: In a sterile package for ease of application.

GENERATES NO HEAT: Contains no heat-generating powder or granules.

FLEXIBLE AND CONFORMABLE:

When applied to the tissue surface.

COLOR-CODED: Simplifies emergency medicine, enabling all responders to select the right product, for the right intervention, at the right time.

SPECIFICATIONS

PN	PDF130
NSN	6510-01-591-7740
Purchase	Direct, Prime Vendor, ECAT, GSA
Packaged dimensions	5.5" x 4.5" x 0.6" (14 x 11.5 x 1.5cm)
Weight	0.7oz (20.5g)
Shelf life*	5.5 years
Misc.	BAA and TAA Compliant; CE marked

* From date of manufacture



visit safeguardmedical.com
call (US) +1 (855) 428-6074 / (Int'l) +44 (0) 1568 613942
email (US) info@safeguardmedical.com / (Int'l) sales.admin@safeguardmedical.com

Devices included in this document may not be cleared or licensed for use or sale in your individual country. Please contact your local distributor for information regarding availability of these products.