



FRESH WHOLE BLOOD, DELIVERED.

Why Fresh Whole Blood?

The transfusion of Fresh Whole Blood (FWB) is the best option for the management of patients in hemorrhagic shock. The ability to rapidly provide warm FWB at or near the point of injury, has proven time and time again to improve outcomes. Whole Blood (WB), and in particular Low-Titer Type-O Whole Blood (LTOWB), is the preferred resuscitation product for the pre-hospital treatment of patients in hemorrhagic shock.¹

Transfusion Family of Products

Our family of FWB transfusion products is designed to **enable, augment, and simplify** the delivery of this life-saving intervention in any environment. Each transfusion set contains everything needed to properly document donor-recipient information and perform the collection-transfusion of one unit of FWB in a compact, sealed package. The Fresh Whole Blood Transfusion Set additionally provides the ability to screen ABO/Rh for the donor and recipient, and collect samples for retrospective testing. Our family of blood transfusion sets are an affordable and sensible solution for training and equipping organizations required to collect/transfuse FWB.



1. Joint Trauma System Clinical Practice Guideline; Whole Blood Transfusion (CPG ID: 21), May 2018, pg. 2.
Devices included in this document may not be cleared or licensed for use or sale in your individual country. Contact your local distributor for availability.
SGM-MKT-FWB-01-002 | © 2026 Safeguard US Operating, LLC | 04/23/2026



TRANSFUSION SETS

DETAILS & SPECIFICATIONS



	FWB Transfusion Set	FWB Donor Set	FWB Recipient Set	LTOWB Transfusion Set	Adv. LTOWB Collection Set	Adv. LTOWB Admin. Set
PN	80-801	80-801-D	80-801-R	80-875	80-820	80-821
NSN	6515-01-657-4750	6515-01-664-0306	6515-01-663-9469	6515-01-684-9074	6515-01-686-2423	6515-01-686-2401
Shelf life	12 months*	12 months*	12 months*	24 months*	24 months*	24 months*
DHA standardized set	✓					
DHA standardized component		✓	✓			
Compatible w/ pre-screened, pre- identified LTOWB donor programs				✓	✓	✓
Exceeds ASBP CPG for walking blood banks	✓	✓	✓			
Includes individually-sealed donor and recipient sets	✓					
Includes blood type and Rh factor screening test		✓	✓		✓	✓
Compatible w/ USMC Valkyrie Program					✓	✓
Includes blood collection tubes for retrospective analysis of donor blood		✓				

FWB = Fresh Whole Blood | LTOWB = Low-Titer Type-O Whole Blood | DHA = Defense Health Agency | ASBP = Armed Services Blood Program | CPG = Clinical Practice Guidelines | *Up to

visit safeguardmedical.com
call +1-855-428-6074
email sales@safeguardmedical.com

SAFEGUARD
MEDICAL™