



WATERJEL®

BURN DRESSINGS & KITS

Be Prepared For Any Burn Emergency

Over thirty years ago, Water-Jel® Technologies introduced water-based gel technology—a revolution in emergency first response for burn injuries. Today, Water-Jel burn products are used by emergency, fire-rescue, industrial safety, and military professionals around the world.

When Seconds Count

The immediate cooling of burns is recommended by EMS protocols and curriculum. Water-Jel burn products **cool burns on contact** while covering and protecting burn injuries. Our dressings consist of a medical-grade non-woven pad, pre-moistened with Water-Jel's **clear, fragrance-free, water-based gel** with hyaluronic acid.

Water-Jel products do not require re-wetting and are suitable for use in the absence of clean, cool running water. Our gel is comprised of 96% water, and formulated to slowly evaporate, ensuring blankets and dressings remain moist and non-adherent to deliver sustained thermal—*not evaporative*—cooling of burns. Safe and easy to use, Water-Jel burn products are an effective option for managing burn injuries when seconds count.

First Response, Safety & Preparedness

Water-Jel Burn Kits are available in a variety of sizes and feature unique product assortments for different work environments:

- Fire Rescue
- EMS
- Manufacturing
- Electric & Gas Utilities
- Foundries
- Law Enforcement
- Military



Specifications

- Dressings available in six sizes:
 - 2" x 6" → 8" x 18"
 - 4" x 4" → 8" x 20"
 - 4" x 16" → 12" x 16"
- Pre-moistened and easy to use—won't stick to burn injuries
- Location-specific burn dressings available for hand and facial burns
- Medical device—no drugs or active ingredients
- Made in the USA

Note: Soft-sided burn kits can be stored on a wall or in a vehicle, behind and under a seat. Hard plastic burn kits and burn stations include mounting holes.

Devices included in this document may not be cleared or licensed for use or sale in your individual country. Contact your local distributor for availability | WJ-PUB-DRE-01 | ©2026 Safeguard US Operating, LLC | 02/27/2026

visit safeguardmedical.com
 call +1-855-428-6074
 email sales@safeguardmedical.com

SAFEGUARD
 MEDICAL™



BURN KIT + BURN WRAP

large, soft-sided

4" x 4" Burn Dressing.....	3EA	3' x 2.5' Fire Blanket	1EA
4" x 16" Burn Dressing...	1EA	4yd x 3" Sterile Gauze....	4EA
8" x 18" Burn Dressing...	1EA	Scissors.....	1PR
12" x 16" Facial Burn Dressing	1EA	Instruction Card	1EA
		Carry Bag.....	1EA

PN.....	EBK2-HA
Count*	3/kits
Dims*	15"x13"x8"
Weight*	22.3lbs
Cube*	0.92ft³



BURN KIT

small, soft-sided

4" x 4" Burn Dressing.....	2EA	4yd x 3" Sterile Gauze....	2EA
4" x 16" Burn Dressing...	2EA	Scissors.....	1PR
12" x 16" Facial Burn Dressing	1EA	Instruction Card	1EA
		Carry Bag.....	1EA

PN.....	EBK1-HA
Count*	6/kits
Dims*	15"x13"x8"
Weight*	13.7lbs
Cube*	0.92ft³



CRITICAL BURN KIT

soft-sided

4" x 4" Burn Dressing.....	2EA	4yd x 3" Sterile Gauze....	2EA
8" x 20" Hand Burn Dressing	1EA	Scissors.....	1PR
12" x 16" Facial Burn Dressing	1EA		

PN.....	BKCK-HA
Count*	6/kits
Dims*	15"x13"x8"
Weight*	19.9lbs
Cube*	0.92ft³



INDUSTRIAL & WELDING BURN KIT

2" x 6" Burn Dressing.....	2EA	8" x 18" Burn Dressing...	1EA
4" x 4" Burn Dressing.....	2EA	Burn-Jel, Unit Dose.....	25EA
4" x 16" Burn Dressing...	1EA	Elastic Bandage.....	1EA
12" x 16" Facial Burn Dressing	1EA	Instruction Card	1EA

PN.....	IWK-HA
Count*	5/kits
Dims*	16"x12"x8"
Weight*	23.3lbs
Cube*	0.89ft³



FOOD SERVICE BURN KIT

2" x 6" Burn Dressing.....	2EA	Burn-Jel, Unit Dose.....	25EA
4" x 4" Burn Dressing.....	2EA	Elastic Bandage.....	1EA
4" x 16" Burn Dressing...	1EA	Instruction Card	1EA
12" x 16" Facial Burn Dressing	1EA		

PN.....	FSK-HA
Count*	5/kits
Dims*	16"x12"x8"
Weight*	12.3lbs
Cube*	0.89ft³



EMERGENCY BURN STATION

large

4" x 4" Burn Dressing.....	2EA	Burn-Jel, Unit Dose.....	10EA
4" x 16" Burn Dressing...	2EA	Burn-Jel, 2oz Bottle.....	1EA
12" x 16" Facial Burn Dressing	1EA	4yd x 3" Sterile Gauze....	2EA
8" x 18" Burn Dressing...	1EA	Scissors.....	1PR

PN.....	EBSL-HA
Count*	5/stations
Dims*	18"x22"x13"
Weight*	12.6lbs
Dimst	17"x21"x2"
Weightt	2.3lbs



EMERGENCY BURN STATION

small

2" x 6" Burn Dressing.....	2EA	Burn-Jel, Unit Dose.....	4EA
4" x 4" Burn Dressing.....	2EA	4yd x 3" Sterile Gauze..	1EA

PN.....	EBSS-HA
Count*	5/stations
Dims*	15.5"x9"x9"
Weight*	2.7lbs
Dimst	15"x8.5"x1.5"
Weightt	0.5lbs



COMPACT BURN KIT

2" x 6" Burn Dressing.....	1EA	Burn-Jel, Unit Dose.....	6EA
4" x 4" Burn Dressing.....	1EA	Scissors.....	1PR

PN.....	BK10-HA
Count*	6/kits
Dims*	16"x12"x8"
Weight*	23.3lbs
Cube*	0.89ft³



COMPACT BURN KIT PLUS+

2" x 6" Burn Dressing.....	1EA	Burn-Jel, Unit Dose.....	6EA
4" x 4" Burn Dressing.....	1EA	4yd x 3" Sterile Gauze....	2EA
4" x 16" Burn Dressing...	1EA	Scissors.....	1PR

PN.....	BK11-HA
Count*	6/kits
Dims*	16"x12"x8"
Weight*	23.3lbs
Cube*	0.89ft³



*Case measurement | †Single-pack measurement

visit safeguardmedical.com
 call +1-855-428-6074
 email sales@safeguardmedical.com

SAFEGUARD
 MEDICAL™